# Compex®

*Fit 5.0*

SP 6.0

SP 8.0

Upute za korištenje

|  |    |
|--|----|
| 1. Objašnjenje simbola                               | 02 |
| 2. Kako elektrostimulacija radi?                     | 03 |
| 3. Kako MI (mišićna inteligencija) tehnologija radi? | 05 |
| 4. Upute   | 07 |
| Sadržaj kompleta i pribor                            | 07 |
| Opis uređaja   | 08 |
| Prvo korištenje                                      | 10 |
| Funkcije uređaja                                     | 11 |
| Punjenje   | 25 |
| 5. Rješavanje problema                               | 29 |
| 6. Održavanje uređaja                                | 33 |
| Garancija  | 33 |
| Održavanje   | 33 |
| Skladištenje/transport i korištenje                  | 33 |
| Odstranjivanje kvarova                               | 33 |
| 7. Tehničke specifikacije                            | 34 |
| Općenite informacije                                 | 34 |
| Stimulacija živaca                                   | 34 |
| RF podaci  | 35 |
| Standardi  | 35 |
| Informacije o elektromagnetskoj stimulaciji          | 35 |
| 8. EMC tablica                                       | 36 |

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## 8. EMC TABLE

The Compex Stimulator requires special precautions regarding EMC and needs to be installed and put into service according to the information provided on the EMC in this manual.

All RF wireless transmission devices can affect the Compex Stimulator. The use of accessories, sensors, and cables other than those specified by the manufacturer, may result in greater emissions or reduce the immunity of the Compex Stimulator.

The Compex Stimulator should not be used beside or stacked on other equipment, if adjacent or stacked use is necessary, one must check the correct functioning of the Compex Stimulator within the context of the configuration used.

### RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC EMISSIONS

The Compex Stimulator is intended for use in the electromagnetic environment specified below. The customer or user of the Compex Stimulator should ensure that it is used in this environment

| EMISSIONS TEST  | COMPLIANCE     | ELECTROMAGNETIC ENVIRONMENT - GUIDE  |
|---|----------------|--|
| RF emissions<br>CISPR 11                                  | Group 1        | The Compex Stimulator uses RF energy only for its internal operation. Consequently, its RF emissions are unlikely to interfere with any adjacent electrical device (radios, computers, telephones etc.). |
| RF emissions<br>CISPR 11                                  | Class B        | Compex Stimulator is suitable for use in any establishment, other than a private dwelling or a place connected directly to the low voltage mains supply which powers residential buildings.              |
| Harmonic emissions<br>IEC 61000-3-2                       | Class A        |  |
| Voltage fluctuations/ emission oscillations IEC 61000-3-3 | Not applicable |  |

## RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC IMMUNITY

Compex Stimulator is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Stimulator must ensure it is used in this recommended environment.

| IMMUNITY TEST  | TEST LEVEL IEC 60601  | OBSERVANCE LEVEL  | ELECTROMAGNETIC ENVIRONMENT - RECOMMENDATIONS   |
|--|---|---|---|
| Electrostatic discharge (DES)<br>CEI 61000-4-2   | ±6 kV at the contact<br><br>±8 kV in air  | ±6 kV at the contact<br>±8 kV in air  | Floors must be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity must be maintained at a minimum of 30%.  |
| Fast transient electrical bursts<br>CEI 61000-4-4  | ±2 kV for power supply lines<br>±1 kV for input/output lines  | ±2 kV for power supply lines<br><br>Not Applicable for input/output lines   | The quality of the electrical power supply should be that of a typical commercial or hospital environment.  |
| Surge<br>CEI 61000-4-5   | ±1 kV differential mode<br><br>N/A  | ±1 kV differential mode<br><br>±2 kV joint mode   | The quality of the power supply should be that of a typical commercial or hospital environment.   |
| Voltage dips, short interruptions and voltage variations on power supply lines<br>CEI 61000-4-11 | <5 % VT (dips >95 % de UT) for 0.5 cycle<br><br><40 % VT (dips >60 % de UT) for 5 cycles<br><br><70 % VT (dips >30 % de UT) for 25 cycles<br><br><5 % VT (dips >95 % de UT) for 5 seconds | <5 % VT (dips >95 % de UT) for 0.5 cycle<br><br><40 % VT (dips >60 % de UT) for 5 cycles<br><br><70 % VT (dips >30 % de UT) for 25 cycles<br><br><5 % VT (dips >95 % de UT) for 5 seconds | The quality of the power supply should be that of a typical commercial or hospital environment. If the Compex Stimulator user requires continuous operation during mains power cuts, it is recommend that the Compex Stimulator is powered by a UPS or a battery. |
| Magnetic field at grid frequency (50/60 Hz)<br>CEI 61000-4-8                                     | 3 A/m   | 3 A/m   | Magnetic fields at the mains frequency should be at a level characteristic of a typical location in a typical commercial or hospital environment.   |

NOTE :VT is the AC supply voltage before application of the test level.

## RECOMMENDATIONS AND DECLARATION BY THE MANUFACTURER CONCERNING ELECTROMAGNETIC IMMUNITY

Compex Stimulator is designed for use in the electromagnetic environment stipulated below. The buyer or user of the Compex Stimulator must ensure it is used in this recommended environment.

| IMMUNITY TEST  | TEST LEVEL IEC 60601   | OBSERVANCE LEVEL           | ELECTROMAGNETIC ENVIRONMENT - RECOMMENDATIONS   |
|--|--|----------------------------|---|
| <p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p> | <p>3 Vrms<br/>150 kHz to 80 MHz</p> <p>3 V/m<br/>80 MHz to 2.5 GHz</p> | <p>3 Vrms</p> <p>3 V/m</p> | <p>Portable and mobile RF communication devices must only be used relative to the Compex Stimulator and its wiring at a distance which is not less than the spacing recommended and calculated using the appropriate equation for the transmitter's frequency.</p> <p>Recommended spacing<br/> <math>d = 1.2 \sqrt{P}</math><br/> <math>d = 1.2 \sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.3 \sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power of the transmitter in watts (W) set by the manufacturer's specifications and where d is the recommended spacing in metres (m).</p> <p>The field intensity of RF fixed transmitters, as determined by an electromagnetic survey a must be less than the observance level to be found in each frequency rangeb. Interference may occur close to any appliance identified by the following symbol:</p>  |

NOTE 1 At 80 MHz and at 800 MHz, the high frequency amplitude is applied

NOTE 2 These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

a The field intensity from fixed transmitters, such as radio telephone base stations (cellular/wireless) and a mobile radio, amateur radios, AM and FM radio transmissions and TV transmissions cannot be predicted with any accuracy. It may therefore be necessary to consider an analysis of the electromagnetic environment of the site to calculate the electromagnetic environment coming from fixed RF transmitters. If the field intensity measured in the environment where the Compex Stimulator is located exceeds the appropriate RF observance level above, the Compex Stimulator should be monitored to ensure it is operating properly. In the event of abnormal operation, new measures may then be imposed, such as realignment or movement of the Compex Stimulator.

b Above the frequency amplitude from 150 kHz to 80 MHz, the field intensity must be < 3 V/m.

## RECOMMENDED SPACING BETWEEN A PORTABLE AND MOBILE COMMUNICATION DEVICE AND THE COMPEX STIMULATOR

The Compex Stimulator is designed for use in an electromagnetic environment in which radiated RF waves are controlled. The buyer or user of the Compex Stimulator can contribute to preventing electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication devices (transmitters) and the Compex Stimulator according to the table of recommendations below and according to the maximum output power of the telecommunication device.

| MAXIMUM TRANSMITTER OUTPUT POWER W | SPACING ACCORDING TO THE FREQUENCY OF THE TRANSMITTER M |  |   |
|------------------------------------|---|--|---|
|                                    | FROM 150 KHZ TO 80 MHZ<br>$D = 1.2 \sqrt{P}$            | FROM 80 KHZ TO 800 MHZ<br>$D = 1.2 \sqrt{P}$ | FROM 800 MHZ TO 2.5 GHZ<br>$D = 2.3 \sqrt{P}$ |
| 0.01                               | 0.12  | 0.12   | 0.23  |
| 0.1                                | 0.38  | 0.38   | 0.73  |
| 1                                  | 1.2   | 1.2  | 2.3   |
| 10                                 | 3.8   | 3.8  | 7.3   |
| 100                                | 12  | 12   | 23  |

In the case of transmitters whose maximum output power is not shown in the table above, the recommended spacing of d in metres (m) can be calculated using the appropriate equation for the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) as set by the transmitter manufacturer

NOTE 1 At 80 MHz and at 800 MHz, the spacing for high frequency amplitude is applied.

NOTE 2 These guidelines may not be appropriate for some situations. Electromagnetic wave propagation is modified by absorption and reflection due to buildings, objects and persons.

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France SAS